



## At a glance: HEPACCS product overview

HEPACCS instrument  
 HEPACCS instrument with HF



Article no.	95-300   95-300-HF
Weight	90 g
Irrigation & aspiration	yes
Scope of delivery	with separate instrument cable (for 95-300) with separate HF instrument cable (for 95-300-HF)
Reprocessing	150 cycles
HF-function (only for 95-300-HF)	International 3-pin socket for connecting an applied part of type CF (according to IEC 60601-1), monopolar coagulation signal with max. 1.4 kVp
Ultrasonic generator	SONOCA 300 FL-version SONOCA 400 FL-version Software 4.72, 5.72, 5.73

HEPACCS sonotrodes



Article no.	617K0014	617K0022	617K0030
Working length	48 mm	41 mm	312 mm
Outer Ø   Inner Ø	2.2 mm   1.4 mm	3.2 mm   2.2 mm	2.6 mm   1.7 mm
Delivery condition	non-sterile	non-sterile	non-sterile
Reprocessing	1 cycle	1 cycle	1 cycle

HEPACCS flues



Article no.	610K0114-10	610K0113-10	617K0300-10
Reprocessing	single-use	single-use	single-use

HEPACCS torque wrench



Article no.	616K0009
Reprocessing	150 cycles

HEPACCS instrument tray



Article no.	616K0200
Reprocessing	reprocessible

HEPACCS combination tubing



Article no.	700S0318-10
Reprocessing	single-use